

The Caledon Label

The technical and safety data included on the Caledon label complements the MSDS, and provides immediate access to essential information.

Product Name

Hazard Symbols

Caledon Item #

First Aid, Precautionary Measures and Risk Phrases

Proper Shipping Name and UN Number

NFPA Diamond

Lot Number

CAS Number

Technical Specifications

Container Size

REACTIF

REAGENT		FW/PM 58.08	
Meets ACS Specifications	Min. 99.5%	Répond aux exigences ACS	
Assay (by GC)		Teneur (par CPG)	
Maximum Allowable	10	Maximum permissible	
Colour (APHA)	0.001%	Couleur (APHA)	
Residue after evaporation	0.001%	Résidu après évaporation	
Solubility in water	Passes test/Conforme	Solubilité dans l'eau	
Titration acid	0.0005 meq/g	Acide titrable	
Titration base	0.002%	Base titrable	
Aldehyde (as HCHO)	0.05%	Aldéhydes (en tant que HCHO)	
Isopropyl alcohol	0.05%	Alcool isopropylique	
Methanol	0.05%	Méthanol	
Substances reducing permanganate	Passes test/Conforme	Substances réductrices du permanganate	
Water (volumetric KF)	0.5%	Eau (KF volumétrique)	

LOT: 57011

CAS No. 67-64-1

4L

Flammable • Inflammable

ACETONE UN 1090

CALEDON LABORATORIES LTD.
Georgetown, Ont., Canada. L7G 4R9
FAX (905) 877-6666